

Product datasheet for SC203037

C21ORF56 (SPATC1L) (NM_032261) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C21ORF56 (SPATC1L) (NM_032261) Human 3' UTR Clone
Symbol:	C21ORF56
Synonyms:	C21orf56
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032261
Insert Size:	249 bp
Insert Sequence:	<p>>SC203037 3'UTR clone of NM_032261</p> <p>The sequence shown below is from the reference sequence of NM_032261. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p>

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GAGGACGGCAAGCCCCTCTTCGCCTGGTGAGCCGCCCCGCGCCCGCCGCTTGCCTGCAGTAAACGCGT
TTGTTCCAACCCGGGGCCGCGGTGCCTCCTGCGCGTCCCCCGGAGGGGAAAGGGCCGCGTCCCCCGCG
CGCGAGGCCAGAGAAGGCCCGCTCCCACCGGTGCTGGGCCCCGACCGCAGCCCGCGCTGCCCGCACC
TGCAGAGTGCTTCTCACCCTCATTAATAATCATCCGTTTGCA
AGCGGACCGACTTACGCGTAAGCGCGCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCC
CAACCTGCCATCACGAGATTTTCGATTCCACCGCCG
  
```

Restriction Sites:	SgfI-RsrII
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_032261.5</u>
Locus ID:	84221


[View online »](#)

MW: 8.6