

## Product datasheet for **SC203069**

### TTC32 (NM\_001008237) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	TTC32 (NM_001008237) Human 3' UTR Clone
Symbol:	TTC32
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001008237
Insert Size:	386 bp
Insert Sequence:	<p>&gt;SC203069 3'UTR clone of NM_001008237 The sequence shown below is from the reference sequence of NM_001008237. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b>=Stop Codon <b>Red</b>=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA<b>GCGATCGCC</b> CAAAGAAGAAATGTTGCAAAAAATTAT<b>TGA</b>TATTTTTAACTTAATGGAAGTATTGATTCATGATCCTTA CATCTGCATCTAGTTATCAGTAATTTAGATATTGAGCTATTTTGATTTATATTTAAGAAATTAATACAT TAGCACTGAAAGTTAAATAGTGTGTTAAAGGTAGTTAATTTTCAGGTTGAATGGGTTTTTTTTAATGAAG TGTAATAATACCAATGTATAAGTGTATATTATTAAATATTATAGTAAAAAGGAATGTGTGGTA TTTTCTTCAGCAAAACTATTTTTGTGATTTTTTTATTCTCAACTTTTTATTTAAAAAATGTTACATCTG CAGAAAAGTTGAAAGTATAATAAACTCTTCAGTTAAATTCA <b>ACGCGT</b>AAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<a href="#">NM_001008237.3</a>



[View online »](#)

Locus ID: 130502

MW: 15.5