

## Product datasheet for **SC203161**

### TSNARE1 (NM\_145003) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	TSNARE1 (NM_145003) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	TSNARE1
ACCN:	NM_145003
Insert Size:	275 bp
Insert Sequence:	>SC203161 3'UTR clone of NM_145003 The sequence shown below is from the reference sequence of NM_145003. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_145003.5</a></u>
Locus ID:	203062
MW:	10

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATCATCATGCCACCTCTGTCCGAAAGTGTGCTACCCGTGCATCGGCCGCTGAACCGTCAGCTGGCTG
TGGTGTGGGCTCTCCAGCCCCGAGCACCCTCCTGGCCCGGAGCTGGGAGCTGATGGAGTCCGCTGCGT
TGTGAGCTGATGCCGTCCTTGGCCTGCAATGGACCCGAGAGCTTACAGCTGCCTGCTGTCTCCCT
GCACTGACCCCACTCGCCCTGTGTGCGCCCACTGAGCACCCGAATAAACCGCTGTGGCTTATCCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```



[View online »](#)