

Product datasheet for **SC203290**

NADE (BEX3) (NM_206915) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: NADE (BEX3) (NM_206915) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: BEX3
Synonyms: Bex; DXS6984E; HGR74; NADE; NGFRAP1
ACCN: NM_206915
Insert Size: 367 bp
Insert Sequence: >SC203290 3'UTR clone of NM_206915
The sequence shown below is from the reference sequence of NM_206915. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CATCATGATGAATTTGCCTTATGCCTTGACTCCTGCCATTTATCATGAGATTAATACTGTGATTCCCG
CTGTTTTCTTTTTCCTTGCATTTTCCTAATATGCCTTTACTGATCCGTTTGTGTGAACCCATGTTAT
TTCCATGTGTCAAGTGGGTCTTGTGTGGCAGCTTCTATTTGAAGATTGCCTTTGCACTCAGTGTAAAGT
TTCTGTCAAGCAGTGTTCACCCATTTGCATGGAAAAATTTAAAGCTAATAAAGCAATTTAAAAAGCAA
TCTATACATTTATTGTCTCATAAAAAATGTATACTTACATATAACAGATGAACTAAAAAAAATCTGAG
ATTAATACTAGAGTGTAA
ACGCGTAAAGCGCCGCGGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_206915.3](#)



[View online »](#)

Summary: May be a signaling adapter molecule involved in p75NTR-mediated apoptosis induced by NGF. Plays a role in zinc-triggered neuronal death (By similarity). May play an important role in the pathogenesis of neurogenetic diseases.[UniProtKB/Swiss-Prot Function]

Locus ID: 27018

MW: 14.5