

Product datasheet for **SC203373**

CDC26 (NM_139286) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CDC26 (NM_139286) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CDC26
Synonyms:	ANAPC12; APC12; C9orf17
ACCN:	NM_139286
Insert Size:	301 bp
Insert Sequence:	>SC203373 3'UTR clone of NM_139286 The sequence shown below is from the reference sequence of NM_139286. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TCATCTCAATTTGGAAGTCTTGAATTTTAGAGATGGATTATCTTGCATGCCAGAGCGCTGGAATGGAAT AAAATGATGGCAGAAGTACAAACCAGATTTAGAGAATTGAGTGCTTGCAAGTCAAGCAGAATGTACCTCC TGCAGAGACAAATCTTCTGCATGAGATTACTGATGCTTCACTTGCAGTCTAAGCTGGAATCCAAACTCT GGTTTGTCTCTTGAATTTGACTCTATAAACTGATCTGATTTTCTGTTTTTAAAAATAAATATATTT TTGAAAAATGTGAGGCATCCTTAA ACGCGT AAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_139286.4</u>



[View online »](#)

Summary: The protein encoded by this gene is highly similar to *Saccharomyces cerevisiae* Cdc26, a component of cell cycle anaphase-promoting complex (APC). APC is composed of a group of highly conserved proteins and functions as a cell cycle-regulated ubiquitin-protein ligase. APC thus is responsible for the cell cycle regulated proteolysis of various proteins. [provided by RefSeq, Jul 2008]

Locus ID: 246184

MW: 11.5