

Product datasheet for **SC203384**

C6ORF154 (LRRC73) (NM_001012974) Human 3' UTR Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | 3' UTR Clones |
| Product Name: | C6ORF154 (LRRC73) (NM_001012974) Human 3' UTR Clone |
| Symbol: | C6ORF154 |
| Synonyms: | C6orf154 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pMirTarget (PS100062) |
| ACCN: | NM_001012974 |
| Insert Size: | 300 bp |
| Insert Sequence: | <p>>SC203384 3'UTR clone of NM_001012974</p> <p>The sequence shown below is from the reference sequence of NM_001012974. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <p>GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC GACAGTCTGTTGGCTGAGACCGAGATGTACTCTCCACTGGGCTCTGCACACCATTTCACTTATCTAT GTCCCGAGCACCTTGCCCGAGATATCAGGGTCAGGCCCTGGGACTTGGGAGGGAAGTGGGGTCAGGGGC TGCATGGGGGCTCTGGCAGCTCCTGGCAGTGTGGTGGGAAGGAAGCTCTGGAAGCTGTGACTGAGCAAC AGCCTTGGGGGGCACTTGAACCCACGGCAATGCCTTTGAAGTGGCAGCTCTGGCTGCAACCCGCTGGC TCGGAAAAGATTTTATGAACTCCA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG</p> |
| Restriction Sites: | Sgfl-MluI |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | <u>NM_001012974.4</u> |


[View online »](#)

Locus ID: 221424

MW: 10.6