

Product datasheet for **SC203418**

Fbx32 (FBXO32) (NM_148177) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Fbx32 (FBXO32) (NM_148177) Human 3' UTR Clone
Symbol:	Fbx32
Synonyms:	Fbx32; MAFbx
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_148177
Insert Size:	2000 bp



[View online »](#)

Insert Sequence:

>SC203418 3'UTR clone of NM_148177

The sequence shown below is from the reference sequence of NM_148177. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGGACTTTATCAACTTGTTCAGTTCGAATCCAGCACATGACAACACTTCAGAAGGGTCCCCCTGC
TGAAGTGGAGAGCTGGAAATATGGCATTGGACACTTCATTTGTAATAGTGTACATTTTAAACATTGGC
TCGAAACTTCAGAGATAAGTCATGGAGAGGACATTGGAGGGGAGAAATGCAGTTGCTGACTGGGAATTT
AAGAATGTGAACCTCTACTAGAATTGGTATGGAAAAGCAAAATACTGTAATAAACTTTTTTCTAAC
AATTTGCCAGCAAGACTATAAGGGCAATAATTCTATTTCCAGCGGTGAAAATGGAGTCTCTTAATGGTC
ACAGAACTCTCTATAGTTCCTAGGAAGAAAAGGCAAACTCAAATACAAAATAGGACGCTTTGTT
TACAATGTGAAAATTTGTTTAGAAAAGAAAAATGAAGAAGGAAAATATATAGGAGAAAATCCTGGGC
TTTGGGGTTTGATCTGGTGTGTTTGTGAGAAGGCGAAGGAAGGCCACCCATCCTAAACCTGCATGGG
CACAGAGCCTTACCCTACAGAGATACCACACAATGGTCTACCTCTAAAAGCGTAGAGTGCCTCTCTGA
CAACGCATGTGATCCAGTAGGCAGTGTCTGTTTTCTCGTAAGTGTTCATCACTATAAAAAGTGATT
TAAAATGGAAAGGGATATTGGAGGCAATCCTCTGGTTTACACTTCCTTTTTCAGGCAAACAGGGCCCTC
CTTAGAGTTATCTGTTCCACCAGAACAACCCAAAAGCACATGTTTCTTTTCTGTGCATGTTGTAGCTC
ATTTGGGGGAAAAGAAGGCCAATTAACAGTCCCAGGGGAGGATTGCTTGAGCCCAGGAGTCTGAGG
TTTGCAGTGAGCTATGATTGTGCCACTGCACTCCAGCCTGTCTCCAAAACAAAACAAAACAAAACAAA
GCCGTCAGAGGAAAGTCTGCTGCAGAAAAGATGATAGAATTGGTCACCCAATCAATAGCTGTGAAGA
CTTCTTATTCCTGCTGAACTTCCTTGGATATCTTCATGGATAAGCCTGTGTCTGTGGACAGAGT
GGGGCCCAAGTGCTGAGGCTCGAAAAGACAGGTTTTACAATCAGCCAATTAACATCTTACCAGTTGCTG
TACCCCTTGCCAAGGGATGAGCTGAGGAGCTTCTTTTCTTAGTCCAACATCAAAATGAGTCCATCT
GAGCTAGGGCTTTTTGCAAAGGTGATTAATAAATCAAGGCTGGAGTTCCAGCCAAATAGGGAGGAGAGG
ACCACAAGTTCATCTTAACTTGCTTCCGGGCTGGGTAGTTAAAACAGGAAGATCCCAAGGGGATCTTA
CGAGGCAAATATCAGCAGGTAAGTGTGGAAGAAAGGAAATCTCAGGACCTAATGAGTTCGTTGTAATAA
ATTCTGGGAATAGATGGGGATCTACTCCAGTTTTTTACTTTTTACAGAATATTGTTGCTTTCCTAC
AACATGTACATATTTGCAGTTGTATATGCTTTTTTTTTTCCAAATAAACTTGTCCACCTGCATGC
CCTTGCAAATAAGTGAAGCAGAAATAGGAACACAGTCCACATTCAGTTGAGGAACAGTGTATCTTTA
AGAGCTGAGCTTTGGGTGACCTGGAAAGGGGAAAGATGGCTAAGCATGGAGAGAAACGAGGCAAGAGA
CAAGCTATGATAACAACCGCTTCAGCCCTGCCCTCAATAGCACACAACCCACATATCAGCTTTCTGA
AGAGAAGGAACCTACTGTTAGTGTCTCCTCACTTTGCAATGTTTGTGCTACGCCAGAATTTCTCCAGTT
TTTTTCATTATCATCCCCCTGAGAAAAAATTACATTGAATTTAAATTTCCCTAATAAGAGAAATTA
ATATGAAAGAATAGGATTTGTTGGTAAGATTGAGCTTTGGAAGGTCACGAACCATTATCTATCTA
ACGCGT AAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
    
```

Restriction Sites:

Sgfl-Mlul

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq:

[NM_148177.3](#)

Summary:

This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbxs class and contains an F-box domain. This protein is highly expressed during muscle atrophy, whereas mice deficient in this gene were found to be resistant to atrophy. This protein is thus a potential drug target for the treatment of muscle atrophy. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jun 2011]

Locus ID:

114907

MW:

76.6