

## Product datasheet for **SC203431**

### TTC9C (NM\_173810) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** TTC9C (NM\_173810) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** TTC9C  
**ACCN:** NM\_173810  
**Insert Size:** 284 bp  
**Insert Sequence:** >SC203431 3'UTR clone of NM\_173810

The sequence shown below is from the reference sequence of NM\_173810. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGCAGCTCTACCTGGGCATGTTTGGTAACAAAGAAGAAAGATGCTCCTCCAGTTGAACTTAGGTGGA
CCATTAACATGCATGAAGGAGAAATCTGAGCCTCAGCAAGAGAAATTAACCCTATACCTCTGACCCAG
GTGGATTTTGTCTAGTTCTGCACAACTTCACTACTTAGACAGTCTGAGTCTTTTTCTGTCTATCC
ATCTGTTTATTTCTATACCTTTCAATACATGTTATTGTTGCAGATATTTGGCTTGAGAAATATAATCAG
AAACATA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** Sgfl-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_173810.4](#)

**Locus ID:** 283237

**MW:** 10.9



[View online »](#)