

Product datasheet for **SC203438**

Dysferlin (DYSF) (NM_001130977) Human 3' UTR Clone

Product data:

| | |
|--------------------|---|
| Product Type: | 3' UTR Clones |
| Product Name: | Dysferlin (DYSF) (NM_001130977) Human 3' UTR Clone |
| Vector: | pMirTarget (PS100062) |
| Symbol: | DYSF |
| Synonyms: | FER1L1; LGMD2B; LGMDR2; MMD1 |
| ACCN: | NM_001130977 |
| Insert Size: | 301 bp |
| Insert Sequence: | >SC203438 3'UTR clone of NM_001130977 The sequence shown below is from the reference sequence of NM_001130977. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GCCATGAAGCTGGTGAAGCCCTTCAGCTGAGGACTCTCCTGCCCTGTAGAAGGGGCCGTGGGGTCCCCT CCAGCATGGGACTGGCCTGCCTCCTCCGCCAGCTCGGGCAGCTCCTCCAGACCTCCTAGGCCTGATTG TCCTGCCAGGGTGGCAGACAGACAGATGGACCGCCCACTCCAGAGTTGCTAACATGGAGCTCTG AGATCACCCCACTTCCATCATTTCCTTCTCCCCAACCCAACGCTTTTTTGGATCAGCTCAGACATATT TCAGTATAAAACAGTTGGAACCACA ACGCGT AAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG |
| Restriction Sites: | Sgfl-MluI |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | <u>NM_001130977.2</u> |



[View online »](#)

Summary:

The protein encoded by this gene belongs to the ferlin family and is a skeletal muscle protein found associated with the sarcolemma. It is involved in muscle contraction and contains C2 domains that play a role in calcium-mediated membrane fusion events, suggesting that it may be involved in membrane regeneration and repair. In addition, the protein encoded by this gene binds caveolin-3, a skeletal muscle membrane protein which is important in the formation of caveolae. Specific mutations in this gene have been shown to cause autosomal recessive limb girdle muscular dystrophy type 2B (LGMD2B) as well as Miyoshi myopathy. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2008]

Locus ID: 8291

MW: 10.9