

Product datasheet for **SC203454**

ALS2CR1 (NIF3L1) (NM_001142355) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ALS2CR1 (NIF3L1) (NM_001142355) Human 3' UTR Clone
Symbol:	ALS2CR1
Synonyms:	ALS2CR1; CALS-7; MDS015
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001142355
Insert Size:	302 bp
Insert Sequence:	>SC203454 3'UTR clone of NM_001142355 The sequence shown below is from the reference sequence of NM_001142355. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC ACTGACAGGGACCCTCTTCAGGTGGTA TAA TTGCAGAAACATCAGGATAACACATTCTACAAATCAGCT GGATGCCAACTTAAATTTGTAACATGAGTCAGTGGGACTGGTGTGCTTCCAGAGAGTGTCTTCGAGGGT ATCATCATTTCCGGTTTGTAACTTATTACCAAATGTTCTATCGCTCGTAAGGTAAAACCTGTAATAT AACTACCATATTAATAACAAATGTTTCATTATAAACTCTAGGAAAGATTGAATAAAATCTGTTTACTTA ACATTCTTAAAAAAAAAAAAAAAAAAAA ACGCGT AAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001142355.1</u>



[View online >](#)

Summary: May function as a transcriptional corepressor through its interaction with COPS2, negatively regulating the expression of genes involved in neuronal differentiation.[UniProtKB/Swiss-Prot Function]

Locus ID: 60491

MW: 11.5