

Product datasheet for SC203579

EBNA1 binding protein 2 (EBNA1BP2) (NM_006824) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	EBNA1 binding protein 2 (EBNA1BP2) (NM_006824) Human 3' UTR Clone
Symbol:	EBNA1 binding protein 2
Synonyms:	EBP2; NOBP; P40
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006824
Insert Size:	289 bp
Insert Sequence:	<p>>SC203579 3'UTR clone of NM_006824</p> <p>The sequence shown below is from the reference sequence of NM_006824. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p>

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AGAGAGAAGATGAAGAACAGAACACACTAATAGCATCTTTGAATACAAAGAACCAAGAAAAAGGAATG
AAGACTCGCAATTTACGACACACTTTGATCCCTTCTGTTGGTGTGATGTTGTAACATTTCTTTCAAT
AAACTAAAGAAAAATTATTAAGGAACACATACCTTTGGTTAAATAGTCTAGACTAAAAGATTGAGAAG
TTACTTTCCATTGCTATCTATTGATAATTTAGACATTGAGTTCAAATTGCCTTCATTTTATGATAAATA
ATGATTTAACTGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_006824.3</u>


[View online »](#)

Summary: Required for the processing of the 27S pre-rRNA.[UniProtKB/Swiss-Prot Function]

Locus ID: 10969

MW: 11.6