

Product datasheet for **SC203598**

HNRPM (HNRNPM) (NM_031203) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	HNRPM (HNRNPM) (NM_031203) Human 3' UTR Clone
Symbol:	HNRPM
Synonyms:	CEAR; hnRNP M; HNRNPM4; HNRPM; HNRPM4; HTGR1; NAGR1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_031203
Insert Size:	290 bp
Insert Sequence:	>SC203598 3'UTR clone of NM_031203 The sequence shown below is from the reference sequence of NM_031203. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC ATTGACGTTTGAATTGATAGAAACGCT TA AGCAGTTGCCTTTTTAAACATCGATACGAGACCTCTGAA TTTGTATTTTTCTTGTTAACCATTTTAAATTTGTTGGCTGGATGTATAAAGATGTTTAAAAAATTCAGT TGCTTTTTGGGGTAATTTGAATTACTTTTTAATGACTGGGGTCCATTTGACTGTTGCATTGAGATT GCAATGTGCGCAATTTTTTTGTAGTTGTGGCATCTTGTGACATCGAATATGACTTTGATAATAAATA CCGGTTCCTGAAAA ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_031203.4</u>



[View online »](#)

Summary:

This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has three repeats of quasi-RRM domains that bind to RNAs. This protein also constitutes a monomer of the N-acetylglucosamine-specific receptor which is postulated to trigger selective recycling of immature GlcNAc-bearing thyroglobulin molecules. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011]

Locus ID:

4670

MW:

11.1