

Product datasheet for **SC203611**

AP2S1 (NM_021575) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	AP2S1 (NM_021575) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	AP2S1
Synonyms:	AP17; CLAPS2; FBH3; FBH0k; HHC3
ACCN:	NM_021575
Insert Size:	310 bp
Insert Sequence:	>SC203611 3'UTR clone of NM_021575 The sequence shown below is from the reference sequence of NM_021575. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CAGCTGCTGATGCTACAGTCCCTGGAGTGAAGGCGAGGCGAGCCCCACCCCGCCCCGCCCTCCTGGA CTCGCCTGCTCGCTTCCCCTTCCAGGCCCGTGGCCAACCCAGCAGTCTTCCCTCAGTGCCTAGGAG GAAGGGACCCAGCTGGGTCTGGGCCACAAGGGAGGAGACTGCACCCCACTGCCTCTGGGCCCTGGCTGT GGCAGAGGCCACCGTGTGTGTCCCGAGTAACCGTGCCGTTGTCTGTGTATGCCATAAGCGTCTGTGCG TGGAGTCCCAATAAACCTGTGGTCTGCCTGGC ACGCGT AAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_021575.5</u>



[View online »](#)

Summary: One of two major clathrin-associated adaptor complexes, AP-2, is a heterotetramer which is associated with the plasma membrane. This complex is composed of two large chains, a medium chain, and a small chain. This gene encodes the small chain of this complex. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]

Locus ID: 1175

MW: 10.4