

Product datasheet for **SC203683**

HAGHL (NM_032304) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: HAGHL (NM_032304) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: HAGHL
ACCN: NM_032304
Insert Size: 289 bp
Insert Sequence: >SC203683 3'UTR clone of NM_032304
The sequence shown below is from the reference sequence of NM_032304. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGGCTCCTGAGTGCAGCCCCACACGACTGAGCCACCCAGACCCTCACAGGGCTGGGGCTGCGTCCCTC
CTCGTGACCTCGGCCAGCTGGACCCACATGAGGGCCACCTCTGGAACCTTCTTCGAGGCCCTGGCCAGC
CATCTGCCAGCCTCGGAGGGTGGGCAACCTGGTGCTTCCCGGGTGGACACACAGGACCCTCAGTGGG
GCCTGTGTGGGCGCCGAGACCTGGGTGTCTGGGAAGTGGGGCACACGGGGCTCCGAACATGAATAAA
GCTTTGAAAGGCC
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_032304.4](#)

Summary: Hydrolase acting on ester bonds.[UniProtKB/Swiss-Prot Function]

Locus ID: 84264

MW: 9.8



[View online >](#)