

Product datasheet for **SC203730**

RPL28 (NM_001136136) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RPL28 (NM_001136136) Human 3' UTR Clone
Symbol:	RPL28
Synonyms:	L28
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001136136
Insert Size:	318 bp
Insert Sequence:	<p>>SC203730 3'UTR clone of NM_001136136 The sequence shown below is from the reference sequence of NM_001136136. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre>GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CGGCCCTTTCCCGGGTCTGGGAGCTGTGATTTTTTACTGTCAGGCAGGAAGAGCGGTAAGTCCATCG CGGCGGGCATCCCTGGCGCCAGGGTGTGGTCTGGGTACCGGCTCCCTCTCGGCCGACTTGTCAGCTC TGTGAGCCGCGCGCTCTGAGCCGTGTCCTCACCTGAAAGTGGAGAAATGAAAAAGGACCTGAACTT CCTCGGTGGTTGTTGAGAGTTAAGGCACGGGGTTGATGTTTTTCAGATGAAATTCTCAAAGCAAGTCAGG GTGGGGATGGATGGTTTCATCCACAGGTGGGAAGATTGAGG ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG</pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001136136.1</u>



[View online »](#)

Summary:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L28E family of ribosomal proteins. It is located in the cytoplasm. Variable expression of this gene in colorectal cancers compared to adjacent normal tissues has been observed, although no correlation between the level of expression and the severity of the disease has been found. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternative splicing results in multiple transcript variants encoding distinct isoforms.[provided by RefSeq, Oct 2008]

Locus ID:

6158

MW:

12.4