

Product datasheet for **SC203743**

Dysadherin (FXVD5) (NM_144779) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: Dysadherin (FXVD5) (NM_144779) Human 3' UTR Clone
Symbol: Dysadherin
Synonyms: DYSAD; HSPC113; IWU1; KCT1; OIT2; PRO6241; RIC
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_144779
Insert Size: 296 bp
Insert Sequence: >SC203743 3'UTR clone of NM_144779
The sequence shown below is from the reference sequence of NM_144779. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG  
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC  
TCCCGTTATGCCGGAATCGTTGCAGGTGAGTCCATCAGAAACAGGAGCTGACAACCTGCTGGGCACCC  
GAAGACCAAGCCCCTGCCAGCTCACCGTGCCAGCCTCTGCATCCCCTCGAAGAGCCTGGCCAGAGA  
GGGAAGACACAGATGATGAAGCTGGAGCCAGGGCTGCCGGTCCGAGTCTCCTACCTCCCCAACCCCTGC  
CCGCCCTGAAGGCTACCTGGCGCCTTGGGGCTGTCCCTCAAGTTATCTCCTCTGCTAAGACAAAAAG  
TAAAGCACTGTGGTCTTTGC  
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA  
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq: [NM_144779.3](#)



[View online >](#)

Summary:

This gene encodes a member of a family of small membrane proteins that share a 35-amino acid signature sequence domain, beginning with the sequence PFXD and containing 7 invariant and 6 highly conserved amino acids. The approved human gene nomenclature for the family is FXD-domain containing ion transport regulator. Mouse FXD5 has been termed RIC (Related to Ion Channel). FXD2, also known as the gamma subunit of the Na,K-ATPase, regulates the properties of that enzyme. FXD1 (phospholemman), FXD2 (gamma), FXD3 (MAT-8), FXD4 (CHIF), and FXD5 (RIC) have been shown to induce channel activity in experimental expression systems. Transmembrane topology has been established for two family members (FXD1 and FXD2), with the N-terminus extracellular and the C-terminus on the cytoplasmic side of the membrane. This gene product, FXD5, is a glycoprotein that functions in the up-regulation of chemokine production, and it is involved in the reduction of cell adhesion via its ability to down-regulate E-cadherin. It also promotes metastasis, and has been linked to a variety of cancers. Alternative splicing results in multiple transcript variants. [RefSeq curation by Kathleen J. Sweadner, Ph.D., sweadner@helix.mgh.harvard.edu., Sep 2009]

Locus ID:

53827

MW:

11.1