

Product datasheet for **SC203796**

RGS10 (NM_001005339) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RGS10 (NM_001005339) Human 3' UTR Clone
Symbol:	RGS10
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001005339
Insert Size:	326 bp
Insert Sequence:	>SC203796 3'UTR clone of NM_001005339 The sequence shown below is from the reference sequence of NM_001005339. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC AAAAGAGCTTCCAGAATTATAACACA TGA GCCCCAAAAAGCCGGGACTGGCAGCTTTAAGAAGCAAA GGAATTTCTCTCAGGACCGTCCGGGTTATCATTGCTTTGTTATTTGTAAGGACTGAAATGTACAAA ACCTTCAATGGGATGTGTGTTTTATTAAGTCTCACCAGTAAATTTGCATGATGGCTAAGCTAACA TAAAAAAGAATAATAA AACTGGAAGT TTTAGTTTACAAAACAGAGATTCCTTCAACTGGACACG TCGAGCATTTTGTAGCTTAATTAACCTCATGTAAGGCCACAAGGTGAAA ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001005339.2</u>



[View online »](#)

Summary:

Regulator of G protein signaling (RGS) family members are regulatory molecules that act as GTPase activating proteins (GAPs) for G alpha subunits of heterotrimeric G proteins. RGS proteins are able to deactivate G protein subunits of the Gi alpha, Go alpha and Gq alpha subtypes. They drive G proteins into their inactive GDP-bound forms. Regulator of G protein signaling 10 belongs to this family. All RGS proteins share a conserved 120-amino acid sequence termed the RGS domain. This protein associates specifically with the activated forms of the two related G-protein subunits, G-alpha_{i3} and G-alpha_z but fails to interact with the structurally and functionally distinct G-alpha subunits. Regulator of G protein signaling 10 protein is localized in the nucleus. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Locus ID:

6001

MW:

12.1