

## Product datasheet for **SC203855**

### Dysferlin (DYSF) (NM\_001130455) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	Dysferlin (DYSF) (NM_001130455) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	DYSF
Synonyms:	FER1L1; LGMD2B; LGMDR2; MMD1
ACCN:	NM_001130455
Insert Size:	301 bp
Insert Sequence:	>SC203855 3'UTR clone of NM_001130455 The sequence shown below is from the reference sequence of NM_001130455. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GCCATGAAGCTGGTGAAGCCCTTCAGCTGAGGACTCTCCTGCCCTGTAGAAGGGGCCGTGGGGTCCCCT CCAGCATGGGACTGGCCTGCCTCCTCCGCCAGCTCGGGCAGCTCCTCCAGACCTCCTAGGCCTGATTG TCCTGCCAGGGTGGCAGACAGACAGATGGACCGCCCACTCCAGAGTTGCTAACATGGAGCTCTG AGATCACCCCACTTCCATCATTTCCTTCTCCCCCAACCAACGCTTTTTTGGATCAGCTCAGACATATT TCAGTATAAAAACAGTTGGAACCACA <b>ACGCGT</b> AAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_001130455.2</a></u>



[View online »](#)

**Summary:**

The protein encoded by this gene belongs to the ferlin family and is a skeletal muscle protein found associated with the sarcolemma. It is involved in muscle contraction and contains C2 domains that play a role in calcium-mediated membrane fusion events, suggesting that it may be involved in membrane regeneration and repair. In addition, the protein encoded by this gene binds caveolin-3, a skeletal muscle membrane protein which is important in the formation of caveolae. Specific mutations in this gene have been shown to cause autosomal recessive limb girdle muscular dystrophy type 2B (LGMD2B) as well as Miyoshi myopathy. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2008]

**Locus ID:** 8291

**MW:** 10.9