

Product datasheet for **SC203868**

CD37 (NM_001774) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CD37 (NM_001774) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CD37
Synonyms:	GP52-40; TSPAN26
ACCN:	NM_001774
Insert Size:	308 bp
Insert Sequence:	>SC203868 3'UTR clone of NM_001774 The sequence shown below is from the reference sequence of NM_001774. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GTCTACAACCGGCTCGCTCGATACCGT TAG GCCCCGCCCTCCCCAAAGTCCCGCCCCGCCCGTCACG TGCGCTGGGCACTCCCTGCTGCCTGTAATATTTGTTAATCCCCAGTTCGCCTGGAGCCCTCCGCCT TCACATCCCCTGGGACCCACGTGGCTGCGTGCCCTGCTGCTCACCTCTCCCACGGGACCTGGGG CTTTCGTCCACAGCTTCTGTCCCATCTGTCCGCTACCACCACCCACAAGATTATTTTCACCCAAA CCTCAAATAAATCCCCTGCGTTTTTGGTAAAA ACGCGT AAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001774.3</u>



[View online »](#)

Summary:

The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell surface glycoprotein that is known to complex with integrins and other transmembrane 4 superfamily proteins. It may play a role in T-cell-B-cell interactions. Alternate splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Locus ID: 951

MW: 11.2