

Product datasheet for **SC203948**

Protein FAM195A (MCRIP2) (NM_138418) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Protein FAM195A (MCRIP2) (NM_138418) Human 3' UTR Clone
Symbol:	Protein FAM195A
Synonyms:	C16orf14; c349E10.1; FAM195A
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_138418
Insert Size:	314 bp
Insert Sequence:	>SC203948 3'UTR clone of NM_138418 The sequence shown below is from the reference sequence of NM_138418. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC TTCCTGGCGAGAATCACCAGCTGTTCCAGTGGCTGCTGGGAGGGGGCGCTGCTACACGGCCGACCTGT CGCCAGGAGAGAAGCATGGCGCCTGCCACCCACTGCGCCTGGGTGGGTGCCGCCACACCTGAAGTG CCAGCATTTGGACTTTTGCACCTTTTTTCCCTTGGCCCGGCTGTCCCAACCAAGCTGCCATGGCCAAG GGCCGAACCCGTCTGACCTCAGCCCTGCTCACTGTGCCAGGGACCAGCGACCAGCCCTGGGGCTGGC AGGGAGGAGCTCCAGGCTAATAAAGTGGAGAACTGTC ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_138418.4</u>



[View online »](#)

Locus ID: 84331

MW: 11.2