

Product datasheet for **SC204053**

Kanadaptin (SLC4A1AP) (NM_018158) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Kanadaptin (SLC4A1AP) (NM_018158) Human 3' UTR Clone
Symbol:	Kanadaptin
Synonyms:	HLC3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018158
Insert Size:	316 bp
Insert Sequence:	>SC204053 3'UTR clone of NM_018158 The sequence shown below is from the reference sequence of NM_018158. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ACCCATCTTAATGACAAGTATGGCTATATGATTGCTTCAGAATCCCAAAGAAAACCTTGTGGACCATGT
GACATGGAATATTTGGGATAATGTATCAAATTGAATGGCCAGAGAAGTTAGATGATTATTTGTAAGAT
CTGGTGACTGGCTTTTCGTTCTGTGTTCTTGGCTTCCTAAATTTATCTGCCCATATGATTCTCATGCAT
TTGATATTTATGTTTAAAGTGTTTATATATGTATGTAAAAAGGGAACCATATGTTTTGAGAATTTGTA
AAGTGAGAGACATGATCCTATTTAAATAAGAAGGCAAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_018158.3</u>


[View online »](#)

Locus ID: 22950

MW: 12.2