

Product datasheet for **SC204144**

TBC1D3 (TBC1D3F) (NM_032258) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: TBC1D3 (TBC1D3F) (NM_032258) Human 3' UTR Clone
Symbol: TBC1D3
Synonyms: TBC1D3; TBC1D3G
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_032258
Insert Size: 345 bp
Insert Sequence: >SC204144 3'UTR clone of NM_032258
The sequence shown below is from the reference sequence of NM_032258. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG  
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC  
GAAAGTTCTCAGTTCCTCCAGGCTTCTAGAAGCATCTGGGCCAGGGCTCATGGCTGGATAATTTCCCT  
AGGCTTAACAACCAAGCAAGCTTCGCATCCTCGTTTTATTTTTGGTTAACTTATGAAAATGTATTAA  
GAAAGAGTGCAGCTCGAGAGAGATTGAGAGATGGAACACACCAGACCCCGATCACAAAGCCAACCATG  
CCCAGCCCCTCCAGCACCCCCAGCCCCACGACCATCGTTTCTGAATTCTGACGACACCGTGAGCCTGCC  
TTTGACTTCAAACCTCATGGAAGGATAACCACCTTCATGTTTTGAAATAAATGTTTCCTGTGAAATGA  
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA  
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-Mlul
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq: [NM_032258.5](#)



[View online »](#)

Summary: Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11 (By similarity).
[UniProtKB/Swiss-Prot Function]

Locus ID: 84218

MW: 13.6