

Product datasheet for **SC204197**

CACNA1G (NM_198382) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CACNA1G (NM_198382) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: CACNA1G
Synonyms: Ca(V)T.1; Cav3.1; NBR13; SCA42; SCA42ND
ACCN: NM_198382
Insert Size: 753 bp
Insert Sequence: >SC204197 3'UTR clone of NM_198382
The sequence shown below is from the reference sequence of NM_198382. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TCCTCTGACCCAGCAGACCTGGACCCCTGAGTCCTGCCCACTTTCCCACTCACCTTTCTCCACTGGGT
GCCAAGTCCTAGCTCCTCCTCTGGGTATATTCCTGACAAAAGTTCCATATAGACACCAAGGAGCGG
AGGCGCTCCTCCCTGCCTCAGTGGCTCTGGGTACCTGCAAGCAGAACTTCAAAGAGAGTTAAAAGCAG
CAGCCCCGGCAACTCTGGCTCCAGGCAGAAAGGAGAGGCCCGGTGCAGCTGAGTTCCCGACACCAGAAG
CTGTTGGGAGAAAAGCAATACGTTTGTGCAAGTCTCTATGTATATTCTATTTTATAAATTAATTGAAT
CTAGTATATGCGGGATGTACGACATTTTGTGACTGAAGAGACTTGTTCCTTCTACTTTTATGTGTCTC
AGAATATTTTGGAGCGAAGGCGTCTGTCTCTTGGCTATTTAACCTAAAATAACAGTCTAGTTATATT
CCCTCTTCTTCAAAGCACAAGCTGGGACCGGAGCACATTGCAGCCCCAACGGTGGCCCATCTTCAGC
GGAGAGCGAGAACCATTTTGGAACTGTAATGTAACCTATTTTTTCTTTAACCTCGTCATATTTTCT
GTAGGGAAAAAAAAAAGAAAAAGAAAAATGAGATTTTACAAGTAAAATGGAACCTTTTTATATATACA
TACATACATATCTATCTATCTATATATATAAAAATAAAGTAAATTTTCTAAATAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_198382.3
Summary:	Voltage-sensitive calcium channels mediate the entry of calcium ions into excitable cells, and are also involved in a variety of calcium-dependent processes, including muscle contraction, hormone or neurotransmitter release, gene expression, cell motility, cell division, and cell death. This gene encodes a T-type, low-voltage activated calcium channel. The T-type channels generate currents that are both transient, owing to fast inactivation, and tiny, owing to small conductance. T-type channels are thought to be involved in pacemaker activity, low-threshold calcium spikes, neuronal oscillations and resonance, and rebound burst firing. Many alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Sep 2011]
Locus ID:	8913
MW:	29.2