

Product datasheet for **SC204219**

CCDC110 (NM_152775) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CCDC110 (NM_152775) Human 3' UTR Clone
Symbol:	CCDC110
Synonyms:	CT52; KM-HN-1; KMHN1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152775
Insert Size:	350 bp
Insert Sequence:	>SC204219 3'UTR clone of NM_152775 The sequence shown below is from the reference sequence of NM_152775. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTAAAGACAGAACTCTCAAGCATCATTAGGATTCCCTCAAGATTGCTAAGAACTCCCTTCGAGAAATA
GCCAATGTTAATTGTAAGTACTGCTAACTGAAGTATGCATATCTAACTCATCTCAAATGACTCATAATGAGGC
GCTATATATTTGCTTTAAAAAATTATGCACAGAATGGTTTGGTGCAGATAACTTAAGCTTTCTCAGCTGG
GAATTTCAGAAGATGGGGGAGAAGTCTACAAGGGAGGGCAGGTACCTGTACATGGAAGGAGCTGAGAT
GTCTAAAGTCACCCACATTTGTCATTTAAAAAAGTATATATAAATAAAAGATACTATTTTGAAGAAGTA
TAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online »](#)

RefSeq: NM_152775.4
Locus ID: 256309
MW: 13.4