

Product datasheet for **SC204221**

MPST (NM_021126) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: MPST (NM_021126) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: MPST
Synonyms: MST; TST2; TUM1
ACCN: NM_021126
Insert Size: 331 bp
Insert Sequence: >SC204221 3'UTR clone of NM_021126
The sequence shown below is from the reference sequence of NM_021126. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATCTCAGAGGGCCGGGGGAAGACCCACTGAGCTGGCAGGACACAGGCGAGCTCAGGTGATGCCGGCC
ACCAGCAATGCCGGCTGGCTGGTAGCTCCGCTTCTGCTTTCACCAAGAGAGTGTTCCTCACTCAACTCAG
GTGGCATTGGGGTGACATCTCAAAGGCCAGGAATCCGTTGACTTGTTGGCTGCCAGTAGGGGCGGGA
GGAAAGGCGGAGGCGAGCCCTGGAGGAGGAGGCCACAACCTCCGAGCTGCCACCTGGTGTGAGCTGG
GGCCCCGCTCCTTTCTGTTTTATTTTTGAGGAAATAAAATAACCAAGTGCTAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_021126.8](#)



[View online »](#)

Summary: This protein encoded by this gene catalyzes the transfer of a sulfur ion from 3-mercaptopyruvate to cyanide or other thiol compounds. It may be involved in cysteine degradation and cyanide detoxification. There is confusion in literature between this protein (mercaptopyruvate sulfurtransferase, MPST), which appears to be cytoplasmic, and thiosulfate sulfurtransferase (rhodanese, TST, GeneID:7263), which is a mitochondrial protein. Deficiency in MPST activity has been implicated in a rare inheritable disorder known as mercaptolactate-cysteine disulfiduria (MCDU). Alternatively spliced transcript variants encoding same or different isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

Locus ID: 4357

MW: 11.6