

Product datasheet for **SC204272**

CLECSF6 (CLEC4A) (NM_194450) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CLECSF6 (CLEC4A) (NM_194450) Human 3' UTR Clone
Symbol:	CLECSF6
Synonyms:	CD367; CLECSF6; DCIR; DDB27; HDCGC13P; hDCIR; LLIR
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_194450
Insert Size:	350 bp
Insert Sequence:	>SC204272 3'UTR clone of NM_194450 The sequence shown below is from the reference sequence of NM_194450. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTTTGTGAGATGATGAAGATCCACTTATGAACTGAACATTCTCCATGAACAGGTGGTTGGATTGGTATC
TGTCATTGTAGGGATAGATAATAAGCTCTTCTTATTTCATGTGTAAGGGAGGTCCATAGAATTTAGGTGG
TCTGTCAACTATTCTACTTATGAGAGAATTGGTCTGTACATTGACTGATTCACTTTTTTCATAAAGTGAG
CATTTATTGAGCATTTTTTCATGTGCCAGAGCCTGTACTGGAGGCCCCCATTTGTGCACACATGGAGAGA
ACATGAGTCTCTCTTAATTTTTATCTGGTTGCTAAAGAATTATTTACCAATAAAATTATATGATGTGGT
GTCTC
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online »](#)

RefSeq: [NM_194450.3](#)

Summary: This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signalling, glycoprotein turnover, and roles in inflammation and immune response. The encoded type 2 transmembrane protein may play a role in inflammatory and immune response. Multiple transcript variants encoding distinct isoforms have been identified for this gene. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region. [provided by RefSeq, Jul 2008]

Locus ID: 50856

MW: 13.3