

## Product datasheet for **SC204345**

### Dermokine (DMKN) (NM\_001035516) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	Dermokine (DMKN) (NM_001035516) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	DMKN
Synonyms:	UNQ729; ZD52F10
ACCN:	NM_001035516
Insert Size:	349 bp
Insert Sequence:	>SC204345 3'UTR clone of NM_001035516 The sequence shown below is from the reference sequence of NM_001035516. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GGCCTGCTGCAGTGGGTGAAGTTTTGGTAGGCAATTTCTTGAACCACCACCGAGGCCCGAAAAGCAC TGGTCGTCAGGGAGCTCCTCCCCTTGGCCCCAGCCTGTGCCAGCCCTGGCCCGGCTGCCACACCTCTG TTTCTAGGCTGGGACCCAGCTTGTCTCTCCTTGTCTTCTCCACTGCACTGTGGTGCTTCAGTGGCC ACCAGCCTCGTCACATACACCAGCATCTTTCTGTACCTCCTCCCTTTGGTGACCTGAAGTCACTGTGAC AGTTCTCCAGGAAGGAGGAGCTTCTACTTTTGTAGTTTCTGTGGAAATAAACATGAATCTTGTTC CCTA ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_001035516.4</a></u>



[View online »](#)

**Summary:** This gene is upregulated in inflammatory diseases, and it was first observed as expressed in the differentiated layers of skin. The most interesting aspect of this gene is the differential use of promoters and terminators to generate isoforms with unique cellular distributions and domain components. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jun 2010]

**Locus ID:** 93099

**MW:** 12.7