

Product datasheet for **SC204352**

Kallikrein 11 (KLK11) (NM_006853) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Kallikrein 11 (KLK11) (NM_006853) Human 3' UTR Clone
Symbol:	Kallikrein 11
Synonyms:	PRSS20; TLSP
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006853
Insert Size:	339 bp
Insert Sequence:	>SC204352 3'UTR clone of NM_006853 The sequence shown below is from the reference sequence of NM_006853. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC TGGATCCAGGAGACGATGAAGAACAAT TAG ACTGGACCCACCCACCACAGCCCATCACCTCCATTCC ACTTGGTGTGGTTCTGTTCACTCTGTTAATAAGAAACCTAAGCCAAGACCCTACGAACATTCT TTGGCCCTCCTGGACTACAGGAGATGCTGTCACTTAATAATCAACCTGGGGTTCGAAATCAGTGAGACC TGGATTCAAATTCTGCCTTGAAATATTGTACTCTGGGAATGACAACACCTGGTTTGTCTCTGTTGTA TCCCCAGCCCCAAGACAGCTCCTGGCCATATATCAAGGTTTCAATAAATATTTGCTAAATGA ACGCGT AAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_006853.3</u>



[View online »](#)

Summary: Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. Alternate splicing and the use of alternate promoters results in multiple transcript variants encoding distinct isoforms which are differentially expressed. [provided by RefSeq, Dec 2016]

Locus ID: 11012

MW: 12.5