

## Product datasheet for **SC204451**

### ZNF438 (NM\_001143771) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF438 (NM_001143771) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ZNF438
Synonyms:	bA330O11.1
ACCN:	NM_001143771
Insert Size:	359 bp
Insert Sequence:	>SC204451 3'UTR clone of NM_001143771 The sequence shown below is from the reference sequence of NM_001143771. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC ATCGAATTTCCAGTGAAGCTGAGAAA <b>TGA</b> GACCCCAAGGCAGCCTGGGGTTAAGGAGAGAGCTCTGCC GCCACCTTCTCAGAGCTTCGTGCTTTATGGTGGTCTTAGTCACAAAGATCAAACAACAGGATTGGT GTGAGTGAACAGAAATGATTTTTGTACATGTTTTATTTTTCTTAACGAAATAAAATATAAGCTCTCGAA GCATATTTTTCTAACTATATGCAGTCTTATTTTAAAATACGTAAGTATTCTTAATATCTTAGACTTTCA GAATAAATAAAAAGATCAAAATTTGAAGATTGAAATTTAGCAGAAATATAGGCTGTTCTTCAAATACAGC TCTCAGTGAACA <b>ACGCGT</b> AAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_001143771.2</a></u>
Summary:	Isoform 1 acts as a transcriptional repressor.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Locus ID: 220929

MW: 13.9