

Product datasheet for **SC204456**

ZNF438 (NM_001143766) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF438 (NM_001143766) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ZNF438
Synonyms:	bA330O11.1
ACCN:	NM_001143766
Insert Size:	359 bp
Insert Sequence:	>SC204456 3'UTR clone of NM_001143766 The sequence shown below is from the reference sequence of NM_001143766. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC ATCGAAGTTCCAGTGAAGCTGAGAAAATGAGACCCCAAGGCAGCCTGGGGTTAAGGAGAGAGCTCTGCC GCCACCTTCCTTCAGAGCTTCGTGCTTTATGGTGGTCTTAGTCACAAAGATCAAACAACAGGATTGGT GTGAGTGAACAGAAATGATTTTTGTACATGTTTTATTTTTCTTAACGAAATAAAATATAAGCTCTCGAA GCATATTTTTCTAACTATATGCAGTCTTATTTTAAAATACGTAAGTATTCTTAATATCTTAGACTTTCA GAATAAATAAAAAGATCAAAATTTGAAGATTGAAATTTAGCAGAAATATAGGCTGTTCTTCAAATACAGC TCTCAGTGAAAACA ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG Restriction Sites: SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001143766.2
Summary:	Isoform 1 acts as a transcriptional repressor.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Locus ID: 220929

MW: 13.9