

Product datasheet for **SC204459**

ZNF438 (NM_182755) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ZNF438 (NM_182755) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ZNF438
Synonyms: bA330O11.1
ACCN: NM_182755
Insert Size: 359 bp
Insert Sequence: >SC204459 3'UTR clone of NM_182755

The sequence shown below is from the reference sequence of NM_182755. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGCGATCGCC
ATCGAACTTTCCAGTGAAGCTGAGAAAATGACCCCAAGGCAGCCTGGGGTTAAGGAGAGAGCTCTGCC
GCCACCTTCCTTCAGAGCTTCGTGCTTTATGGTGGTCTTAGTCACAAAGATCAAACAACAGGATTGGT
GTGAGTGAACAGAAATGATTTTTGTACATGGTTTTATTTTCTTAACGAAATAAAATATAAGCTCTCGAA
GCATATTTTTCTAACTATATGCAGTCTTATTTTAAAATACGTAAGTATTCTTAATATCTTAGACTTTCA
GAATAAATAAAAAGATCAAAATTTGAAGATTGAAATTTAGCAGAAATATAGGCTGTTCTTCAAATACAGC
TCTCAGTGAAAACA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_182755.3](#)

Summary: Isoform 1 acts as a transcriptional repressor.[UniProtKB/Swiss-Prot Function]



[View online >](#)

Locus ID: 220929

MW: 13.9