

Product datasheet for **SC204478**

ZNF195 (NM_001130519) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ZNF195 (NM_001130519) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ZNF195
Synonyms: HRF1; ZNFP104
ACCN: NM_001130519
Insert Size: 1221 bp
Insert Sequence: >SC204478 3'UTR clone of NM_001130519
The sequence shown below is from the reference sequence of NM_001130519. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTGACTGTACATGAAAGCATTTCATACTTGAAGAAAAATAAACAATATAAAAAATAGGCAAAGCCGTTA
GTATCTGCTCGCATCCCACTTTACATCAGAGTTCAGACTTAATAAAGTTATCATAAATGTAATTACTAT
TGAAAGACCTTTTCATGAAATATAAGTCTCCAGAGCACACAAGAGTATTTCTTCTGAAAAAATGTTACA
ATTATACTAGATGGAAGACTTGCAATGCTTAACTTTTCAAAAAATTTGTAGAGAAACCCAACAAA
TCTCATAAATGTGGAATAACATTTGTTCAAAAATGATAGCTTAGAAAACACTAGAGTTCATACTAAAAG
ATATTTTGCAAATACACTAAATGGGAAAAATAAGCAAGGACACGTTAGAAAATGACACTGAATCCCTGT
GTTGAGTCTACAGAGACATCTAGACTTAAAATAGATCATTATATTTAAGTTTAAAAGGAGGATGAGAC
TGAGTTTTGTAGAGTTATAATTACATTCAAAATATACTTTTTTGCATTGAAAAAATGTTAAATTTTG
AAAAGCGAATATTGATGTCATTCTATTCTCAAATCAGTAGATGCTGTGTCTTCTATTCTAGTGCTGATG
TGAAAACACATGGTCAGTTGTTGCTGCGTCAGAGATATGAGAAATCTTTTCTCTTAGGTGTGCATCAC
TGATATACTTGTCCGTGAAAGGTGAAAGACACTGAAATGTAAGATGCGTGAGTAAATCTTGGTAGAGA
GGCTCTTTGTGGTTGTTTTATATACTTGTAAAGTATTTATGAGATAGGTGTTTAGAATAACTCTTCT
GTTAGGTCACCTGATCAGACCGTTCCCTTCCCTCCCAAGGACTTACAATACGGTCCCTTGTGCTC
TCTGCACAGCTACCCAGGGCAAAAAACAACCCCTTCACTGATCCCTCCAGTAACTGTGCAGACAG
TTACAGGATGCGGTTAACATGCTGTTCACCTCGCATAACAAAGCTGGCAAAAAACATCTCCAGGATGC
AGACAAGACACCTGCACCCTCGACTAAGCTCCCCACCCCAACCCAGTTCTCCTGCACCCCAACTCAG
CTCCCGCACCCCGACCTAGTTCTGGCCCTATAAAAACTGCTATAGTCTGTAAGCAGGGCTGCCTCCTC
TAACGTGGTGGAGCAGCAAGCAGCTCAATAAAGCTTGGCTTGCCTGA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAAATGACCGACCAAGCAGCAGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001130519.3
Summary:	This gene encodes a protein belonging to the Krueppel C2H2-type zinc-finger protein family. These family members are transcription factors that are implicated in a variety of cellular processes. This gene is located near the centromeric border of chromosome 11p15.5, next to an imprinted domain that is associated with maternal-specific loss of heterozygosity in Wilms' tumors. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2012]
Locus ID:	7748
MW:	46.4