

Product datasheet for **SC204499**

PORCN (NM_203475) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PORCN (NM_203475) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PORCN
Synonyms: DHOF; FODH; MG61; PORC; PPN
ACCN: NM_203475
Insert Size: 362 bp
Insert Sequence: >SC204499 3'UTR clone of NM_203475

The sequence shown below is from the reference sequence of NM_203475. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TGCTGGATCTTCTACCGTCTCATAGGCTGAGGCACATCTGTGGACCCTCATAACCCTCTTAAGACCCT
CTCAGGGTGCCACTGATGGGGGATGAGGGAAGGCCCTCTCTACTCCTTGACCCCTCCATCCTTGAC
CCCAACACCTCAACACACACACACACACACACACAAAATCACACCATTTTCATGCCTGTCAATCCC
CACCCACCAAGGCAGGAAGGGGGTGGTGCCTGCTGGGGCCTAGAGGGGGATGCTTGGGGAAACAGA
GAAGGGGAGATCCAGGGCCTCCCGCCTTCTTCTTTTATATAAATTTGTTATTGTCAAATAA
AAGTAGGAAATATTCAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_203475.3](#)



[View online »](#)

Summary: This gene belongs to the evolutionarily conserved porcupine (Porc) gene family. Genes of the porcupine family encode endoplasmic reticulum proteins with multiple transmembrane domains. Porcupine proteins are involved in the processing of Wnt (wingless and int homologue) proteins. Disruption of this gene is associated with focal dermal hypoplasia, and the encoded protein has been implicated in cancer. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Aug 2013]

Locus ID: 64840

MW: 13.5