

Product datasheet for **SC204544**

PSME1 (NM_176783) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PSME1 (NM_176783) Human 3' UTR Clone
Symbol: PSME1
Synonyms: HEL-S-129m; IFI5111; PA28A; PA28alpha; REGalpha
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_176783
Insert Size: 348 bp
Insert Sequence: >SC204544 3'UTR clone of NM_176783
 The sequence shown below is from the reference sequence of NM_176783. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CACTCCCTGACCCTGCAGGCTAGGGGTTAAGGGTGACAAAGCTCAGCTTCTCCACAAGGCTAGAAATGG
GGCACAGAGCCACTGGAGGCCTCTGACTGACCTCTACTCCCTGGCCCTGTAGGCTGTGTTATATGACAT
CATCCTGAAGAATTTCGAGAAGCTCAAGAAGCCCAGGGGAGAAACAAAGGAATGATCTATTGAGAGCC
CTCTCTCCATTCTGTGATGAGTACAGCAGAGACCTTCTGCTTTTTACTGGGGACTCCAGATTTTCCC
CAAACCTGCTTCTGTTGAGATTTTCCCTCACCTTGCTCTCAGGCACAATAAATATAGTTATACCACT
GCC
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online »](#)

RefSeq: [NM_176783.3](#)

Summary: The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core and a 19S regulator. The 20S core is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An essential function of a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. The immunoproteasome contains an alternate regulator, referred to as the 11S regulator or PA28, that replaces the 19S regulator. Three subunits (alpha, beta and gamma) of the 11S regulator have been identified. This gene encodes the alpha subunit of the 11S regulator, one of the two 11S subunits that is induced by gamma-interferon. Three alpha and three beta subunits combine to form a heterohexameric ring. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]

Locus ID: 5720

MW: 12.7