

Product datasheet for **SC204597**

CD72 (NM_001782) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CD72 (NM_001782) Human 3' UTR Clone
Symbol:	CD72
Synonyms:	CD72b; LYB2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001782
Insert Size:	371 bp
Insert Sequence:	>SC204597 3'UTR clone of NM_001782 The sequence shown below is from the reference sequence of NM_001782. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC GAGATGACAGCTTTCAGGTTTCCAGAT TAG GACAGTCCTTGCCTGAGTTGACACTCATGCCAACAAG AACCTGTGCCCTCCTTCTAACCTGAGGCCTGGGTTTCTCAGACCATCTCCTTATTCTGGGCAGTG CCCAGCCACCGGCTGACCCACACCTGACACTTCCAGCCAGTCTGCTGCCTGCTCCCTCTTCTGAAACT GGACTGTTCTGGGAAAAGGGTGAAGCCACCTCTAGAAGGGACTTTGGCCTCCCCCAAGAACTTCCCA TGGTAGAATGGGGTGGGGGAGGAGGGCGCACGGGCTGAGCGGATAGGGGGCGCCCGGAGCCAGCCAGGC AGTTTTATTGAAATCTTTTTAAATAA ACGCGT AAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online »](#)

RefSeq: [NM_001782.3](#)
Summary: Plays a role in B-cell proliferation and differentiation.[UniProtKB/Swiss-Prot Function]
Locus ID: 971
MW: 13.6