

Product datasheet for **SC204680**

SLC44A3 (NM_152369) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC44A3 (NM_152369) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	SLC44A3
Synonyms:	CTL3
ACCN:	NM_152369
Insert Size:	354 bp
Insert Sequence:	>SC204680 3'UTR clone of NM_152369

The sequence shown below is from the reference sequence of NM_152369. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGAACAGAACTCCAGGCCATTGTGAGATAGATACCCATTTAGGTATCTGTACCTGAAAAACATTTTCCTT
CTAAGAGCCATTTACAGAATAGAAGATGAGACCACTAGAGAAAAAGTTAGTGAATTTTTTTTTAAAGAC
CTAATAAACCCCTATTCTTCCTCATTGTCTTTGTCATTATTGTTTGACCAGGTAACAATACTGGAACAT
ATTAGTTTACCTTTTTTTGTACAAATTAGGACAGAAAACTCTTCTAAAACCATGTTTATATGCATCAA
CTTACAAAGTACACTATGTAAGAACTGAGGTAAGTTTGTAAAGTGCACAATAATAAACCTTTTTTA
AGATAAGGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_152369.5</u>
Locus ID:	126969



[View online >](#)

MW: 14.1