

Product datasheet for **SC204736**

UBAP2L (NM_001127320) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	UBAP2L (NM_001127320) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	UBAP2L
Synonyms:	NICE-4; NICE4
ACCN:	NM_001127320
Insert Size:	382 bp
Insert Sequence:	>SC204736 3'UTR clone of NM_001127320 The sequence shown below is from the reference sequence of NM_001127320. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TACAAGCATTCTGGACGGCTGAGAGCTAATTTGGCCCAAGGCTGGGGGCTGTGTTTGTGTGTGTGTA TAAATTTGCACTGAAGTCTTGTTCAGAAACCAGACCACTGAGGAGAGCCTGCTGAGCTGAGGCCATGG CCTGCGTGGCTTGGGAAATGAGTTGGTGATACCTTCTGGGCTTTGAACTTGCCCTCCCCATTTTC CCTCTCCCCATGTGTCTGACCCTGTCTTACCCATTTCAAGTTCAAGCGGTGCAGCACCTTCGAAGCAT CAATGCACACACCTGCTGTTGCTTTTGTATTCTGGAAGGCATGTAGTTTCAACTTGTAAACAAAATATT TGTAGTCTTCAATAAACTGTGGTATTTCTTTAGCTAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG Restriction Sites: SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001127320.3</u>



[View online »](#)

Summary: Plays an important role in the activity of long-term repopulating hematopoietic stem cells (LT-HSCs).[UniProtKB/Swiss-Prot Function]

Locus ID: 9898

MW: 14.1