

Product datasheet for **SC204800**

NPD014 (RSRP1) (NM_020317) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NPD014 (RSRP1) (NM_020317) Human 3' UTR Clone
Symbol:	NPD014
Synonyms:	C1orf63; NPD014
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020317
Insert Size:	370 bp
Insert Sequence:	>SC204800 3'UTR clone of NM_020317 The sequence shown below is from the reference sequence of NM_020317. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AGTCCATATGGACTGTGGATACCTATCAAAAGAAGAAACTGATGGCTAAGTTTGCATGAAAAGTCA
CTTTATTGCAAGTTAGTGTCTAGCATTATCCATCCCTTTGAGCCATTCAGGGGTACTTGTGCATTT
AAAAACCAACACAAAAAGATGTAATACTTAACACTCAAATATTAACATTTTAGGTTTCTCTGCAGAT
ATGAGAGATAGCACAGATGGACCAAAGTTATGCACAGGTGGGAGTCTTTGTATATAGTTGTAATAT
TGTCTTGTTATGTAAAAATGAAATTTTTAGACACAGTAATTGAACTGTATTCCTGTTTGTATATTT
AATAAATTTCTTGTTTTCATTCTTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_020317.5
Locus ID: 57035
MW: 14.5