

Product datasheet for SC204822

IDNK (NM_001001551) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	IDNK (NM_001001551) Human 3' UTR Clone
Symbol:	IDNK
Synonyms:	bA522I20.2; C9orf103; hGntK
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001001551
Insert Size:	373 bp
Insert Sequence:	>SC204822 3'UTR clone of NM_001001551 The sequence shown below is from the reference sequence of NM_001001551. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC ACAATTATGGAACCCCTAAAAATGAAATGCAATGATTTTGTATCAGTGGTCCAAACAGAACTAAGCAT AAATCATTGTGCCATCCCAAACCTCGTTCCAGCCGCTTGCCCACTAGATTCTAAATGTTTCTAAAG GCAAAACCCAATGTGTCAAGACAGACTTGTAGGTGTAATTTAGGAATTATGCTGGTTCATCAGGAA GCAGAGGGGGAGTTTTAAAGTCAAGCTTAAATTGAAGTTTAAATTCATCTATAACCAAATCAAATGAT CAGAGGAAATTCTGTAATCAATGCTGGAAATCGTTACATTGTTTAGAACATTCTTGCTCATGCCTGTAT TTGCACAAATAAATGAACTTCGCTGTC ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_001001551.4

Locus ID: 414328

MW: 14.3