

Product datasheet for **SC204836**

SLC38A10 (NM_138570) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC38A10 (NM_138570) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	SLC38A10
Synonyms:	PP1744
ACCN:	NM_138570
Insert Size:	375 bp
Insert Sequence:	>SC204836 3'UTR clone of NM_138570

The sequence shown below is from the reference sequence of NM_138570. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TACCGTTTTATCTCCCTTGCTCCCTTTAAACACCTGCGGCATTGTCACTTGGCTTGTGATTACCATT
TCGTCTTTAAAGCAGTTTGTAAACAGTTTAACTCGGAGATTTAAAGCCAGCCAGGTGTTGGCGT
GCCTGCACCCAGATGTGGAGTTTATTTAAAGGCTGTGCCCAAGGGGACCCAGTTTGGGGTGGCTGT
GCCCGCCCTGGTCTAGCATAGCCCTGGCCAGAGCACTGCTCCACACCCTCGGCCTCGGAGCTCGG
CTGTGTTTCATCTACAGGATCTGTACCTGTCCCTCCAGGAAAGCTCCTTACGCCTTTGCTCTCCTGAG
AATAAACTTGCTCTGAATCCTCCTCCAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_138570.4
Summary:	Putative sodium-dependent amino acid/proton antiporter.[UniProtKB/Swiss-Prot Function]



[View online >](#)

Locus ID: 124565

MW: 13.5