

Product datasheet for **SC204843**

EHBP1L1 (NM_001099409) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	EHBP1L1 (NM_001099409) Human 3' UTR Clone
Symbol:	EHBP1L1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001099409
Insert Size:	378 bp
Insert Sequence:	>SC204843 3'UTR clone of NM_001099409 The sequence shown below is from the reference sequence of NM_001099409. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA CGGATCGCC AGCCGGCGGGAGCGCTGCGTGCTGAGCT TGA GGCCGCCGGCCCGGGTGGCCATAACTTCTCGCGTCCCC GGCGTCCGCCCGCCCGCCCGGGCCTGCGCTGCGGACGACCCGGCCGTCCCGGAGGCCGCGCGGTGTCCG CTAGGGGCCCGCCGCGCCCTTCCCGTACAGGGCAGGGCGGATCCCCGACCCACGGGCGGGCGGCCG CCGTATTTATTTGTCACCGAGGGTGTGTGCGCGCTCGCGGGGTGCGGGGTCTCCCGACGGCACGG CCGGGCCGCGCGCTCGGGGAGAGGGATGCTGGGCACTACCGCCCGCGCTGGCTTGCCTCTCTGTTT TCCAGAGCAATAAAGTTGGACGAGACTACCCCA ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001099409.3</u>



[View online »](#)

Summary: May act as Rab effector protein and play a role in vesicle trafficking.[UniProtKB/Swiss-Prot Function]

Locus ID: 254102

MW: 13.3