

## Product datasheet for **SC204848**

### **DCTN1 (NM\_004082) Human 3' UTR Clone**

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** DCTN1 (NM\_004082) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** DCTN1  
**Synonyms:** DAP-150; DP-150; P135  
**ACCN:** NM\_004082  
**Insert Size:** 374 bp  
**Insert Sequence:** >SC204848 3'UTR clone of NM\_004082

The sequence shown below is from the reference sequence of NM\_004082. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CACCAGCTTCACAGTCGCTCATCTCCTAAGCACTCCTTTCCCTGTGTCCCTTCGACCCTCAGCCC
TCTGGTGCCGCTCTGCCGATGCACAGCCACCTCAGCCAGCCCCAGGTAGAAACGTGGGTTAAGCTCT
TCTGCCCCGTTAGCTTCACTCCACCCTTTACGCGTCCTGCCCTTACCTTGACCCGGGTTCCCCC
ACTCCCATTCCCTGGCCTCTGCCATAATTTGTTGTTCAACTGCTCCCTCCTTCTGAGGGGCTCAGGG
CTTGTGGGGGTAGGCTGAGACCCACCACAAAGGTTAAGTGAGGTCCCCTTGATTGAGGACTTCACC
CCTTGATTAAGCAACTTCTGCTTCAGTG
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_004082.5](#)



[View online »](#)

<b>Summary:</b>	This gene encodes the largest subunit of dynactin, a macromolecular complex consisting of 10 subunits ranging in size from 22 to 150 kD. Dynactin binds to both microtubules and cytoplasmic dynein. Dynactin is involved in a diverse array of cellular functions, including ER-to-Golgi transport, the centripetal movement of lysosomes and endosomes, spindle formation, chromosome movement, nuclear positioning, and axonogenesis. This subunit interacts with dynein intermediate chain by its domains directly binding to dynein and binds to microtubules via a highly conserved glycine-rich cytoskeleton-associated protein (CAP-Gly) domain in its N-terminus. Alternative splicing of this gene results in multiple transcript variants encoding distinct isoforms. Mutations in this gene cause distal hereditary motor neuropathy type VIIB (HMN7B) which is also known as distal spinal and bulbar muscular atrophy (dSBMA). [provided by RefSeq, Oct 2008]
<b>Locus ID:</b>	1639
<b>MW:</b>	13.4