

## Product datasheet for **SC205112**

### **RNF113B (NM\_178861) Human 3' UTR Clone**

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** RNF113B (NM\_178861) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** RNF113B  
**Synonyms:** bA10G5.1; RNF161; ZNF183L1  
**ACCN:** NM\_178861  
**Insert Size:** 408 bp  
**Insert Sequence:** >SC205112 3'UTR clone of NM\_178861

The sequence shown below is from the reference sequence of NM\_178861. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTTCAGGCTGCAGAAGGAAAAAGAGATAAGAAAGACTCAAGTCTGTATGTACACTGGGAGAAAGCAGC
CTGTTGAGTAAGATGAGAGAGATGACAGTCTTCTGCAATCCCATGCCTGAGTATGTAACAAAAGAAT
TGAATCCGGATCTCAAAGAGGTATGTGCACTCCCATGTTCACTGCAGAATTATTGCAAAAGCCAAGA
GGTGTAACAACCCAAGTGTCCATGGATGGATGAATGGACAAACAAAATGCGATACACACACAATGG
AATGCTATTAAGCCTTAACGGGAAGGAAATCCTTTTAAAGATTTAACTGATCCTTTTAAAGGATCACAT
GCTACAAC TAGAATGAACCTTGAAGACGACACTGCTAAGTAAAATAAACTAGTCACAAAATGA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_178861.5](#)

**Locus ID:** 140432



[View online >](#)

MW: 15.2