

Product datasheet for SC205137

Tartrate Resistant Acid Phosphatase (ACP5) (NM_001111034) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Tartrate Resistant Acid Phosphatase (ACP5) (NM_001111034) Human 3' UTR Clone
Symbol:	Tartrate Resistant Acid Phosphatase
Synonyms:	HPAP; TRACP5a; TRACP5b; TRAP; TrATPase
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001111034
Insert Size:	380 bp
Insert Sequence:	<p>>SC205137 3'UTR clone of NM_001111034</p> <p>The sequence shown below is from the reference sequence of NM_001111034. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC ACCAGGCTGCCGAGGCGAGCCAGGCCCTGAACTCCCATGACTGCCAGCTCTGAGGCCCGATCTCCACT GTTGGGTGGTGGGCCCTGCCGGGACCCTGCTCACAGGCAGGCTTTTCTCCAACCTGTGGCGCTGCAG CAGGGCAGGAAGGGGAACACAGCTGATGAAGTGGTGCCACATGACCCTGTGGCACAGATGCCAC GTATGTGAAACACACATGGACATGTGTCCAGCCACAGTGTATGCTCTGTGGCTGGCTCACCTTTGCT GAGTTCCGGGTGCAATGGGGGAGGGAGGGAGGAAAGCTTCTCCTAAATCAAGCATCTTCTGTTAC TGATGTTCAATAAAGAATAGTTGCCAAGGCTGAA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: [NM_001111034.3](#)

Summary: This gene encodes an iron containing glycoprotein which catalyzes the conversion of orthophosphoric monoester to alcohol and orthophosphate. It is the most basic of the acid phosphatases and is the only form not inhibited by L(+)-tartrate. [provided by RefSeq, Aug 2008]

Locus ID: 54

MW: 13.8