

Product datasheet for **SC205215**

HSBP1L1 (NM_001136180) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	HSBP1L1 (NM_001136180) Human 3' UTR Clone
Symbol:	HSBP1L1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001136180
Insert Size:	401 bp
Insert Sequence:	>SC205215 3'UTR clone of NM_001136180 The sequence shown below is from the reference sequence of NM_001136180. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001136180.2</u>



[View online »](#)

Locus ID: 440498

MW: 15.4