

Product datasheet for **SC205219**

PARP3 (NM_005485) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PARP3 (NM_005485) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PARP3
Synonyms: ADPRT3; ADPRTL2; ADPRTL3; ARTD3; IRT1; PADPRT-3
ACCN: NM_005485
Insert Size: 417 bp
Insert Sequence: >SC205219 3'UTR clone of NM_005485
The sequence shown below is from the reference sequence of NM_005485. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTGCGCTACCTGCTGGAGGTCCACCTCTGAGTGCCCGCCCTGTCCCCGGGGTCTGCAAGGCTGGACT
GTGATCTTCAATCATCTGCCATCTCTGGTACCCTATATCACTCCTTTTTTTCAAGAATACAATACG
TTGTTGTTAACTATAGTCACCATGCTGTACAAGATCCCTGAACTTATGCCTCCTAACTGAAATTTGTA
TTCTTTGACACATCTGCCAGTCCCTCTCTCCAGCCCATGGTAACCAGCATTTGACTCTTTACTTGT
ATAAGGGCAGCTTTTATAGTTCCACATGTAAGTGAGATCATGCAGTGTTTGTCTTTCTGTGCCTGGCT
TATTTCACTCAGCATAATGTGCACCGGGTTCACCCATGTTTTATAAATGACAAGATTTCTCCTTTTT
TAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_005485.6](#)



[View online »](#)

Summary: The protein encoded by this gene belongs to the PARP family. These enzymes modify nuclear proteins by poly-ADP-ribosylation, which is required for DNA repair, regulation of apoptosis, and maintenance of genomic stability. This gene encodes the poly(ADP-ribosyl)transferase 3, which is preferentially localized to the daughter centriole throughout the cell cycle. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Locus ID: 10039

MW: 16