

Product datasheet for SC205223

RP1P8 (RUNDC3A) (NM_001144826) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RP1P8 (RUNDC3A) (NM_001144826) Human 3' UTR Clone
Symbol:	RP1P8
Synonyms:	RAP21P; RP1P-8; RP1P8
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001144826
Insert Size:	387 bp
Insert Sequence:	>SC205223 3'UTR clone of NM_001144826 The sequence shown below is from the reference sequence of NM_001144826. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC GGTCCCATCGGAAGCTCAGAGCCAAATATGTGGCTCCCTTCGAGCGAATGCCCAGGACTTCAACGCATG CACTTTGTGTTGACCTCATCCCTGGCTTACCTTGGTTTTTCCCATCCTAGTTCTCCCTATGCCTGAAT ATCCTGTCTTTTCTTTTATAAGCAACCACACTGTATTGGATGACCCTAGATCTTCTTTGAGACAAGG CAGGCTGTGGCCATGTAGCCCCATCACACTGTGTTTGTGATTGTCTGTGTCTGTCTCCCCACCAGA CTGTGAGCTCCATGAGGGCAGGGACCGTGTCTTGTCCGTTCTCTGTATCCCCAGTGCTTGAACAGAGC GAGTGCTCACTGTGTATTTAATAAATGGACAAAGAGAGAGGA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_001144826.2

Summary: May act as an effector of RAP2A in neuronal cells.[UniProtKB/Swiss-Prot Function]

Locus ID: 10900

MW: 14.4