

Product datasheet for SC205224

RPIP8 (RUNDC3A) (NM_006695) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RPIP8 (RUNDC3A) (NM_006695) Human 3' UTR Clone
Symbol:	RPIP8
Synonyms:	RAP2IP; RPIP-8; RPIP8
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006695
Insert Size:	387 bp
Insert Sequence:	>SC205224 3'UTR clone of NM_006695 The sequence shown below is from the reference sequence of NM_006695. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GGTCCCATCGGAAGCTCAGAGCCAAATAGTGGCTCCCTTCGAGCGAATGCCCAGGACTTCAACGCATG
CACTTTGTGTTGACCTCATCCCTGGCTTACCTTGGTTTTTCCCATCCTAGTTCTCCCTATGCCTGAAT
ATCCTGTCTTTTCTTTTATAAGCAACCACACTGTATTGGATGACCCTAGATCTTCTTTGAGACAAGG
CAGGCTGTGGCCATGTAGCCCCATCACACTGTGTTTGTGATTGTCTGTGTCTGTCTCCCCACCAGA
CTGTGAGCTCCATGAGGGCAGGGACCGTGTCTTGTCCGTTCTCTGTATCCCCAGTGCTTGAACAGAGC
GAGTGCTCACTGTGTATTAAATGACAAAGAGAGAGGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: [NM_006695.5](#)

Summary: May act as an effector of RAP2A in neuronal cells.[UniProtKB/Swiss-Prot Function]

Locus ID: 10900

MW: 14.4