

Product datasheet for **SC205238**

VWCE (NM_152718) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: VWCE

Synonyms: URG11; VWC1

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_152718

Insert Size: 419 bp

Insert Sequence: >SC205238 3'UTR clone of NM_152718
The sequence shown below is from the reference sequence of NM_152718. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GCTTCTCGGGGGGAAGAGTCCACCATGTAAAGGAGTCACTGTGTCCGGGAGACTCTGGAGAGAGGACCT
CTGCCAGTGGCCAGGGTGTGTGCAGGGCAGCTCCAAGGATGAACCTGGTGGGGATGCCTGGGCTCCCT
CCTGCAGGGGCCCTGGTGAGGATGGAAGACCCCAAGGCTGGATGTAACCTTGTCCCAAGAAGTGTTT
GGAATGTGCTGTAAGAATGGAGGAAGTCGTTTCCACTGTCAGCATCCTCCCTGGACCGCGTGGCTGGCT
CATCTTTTGAGAAGGGTTGGGACTGCCAAGTTCTCCTGGAGGAAGAGTTGCGTCCGGCTGGGATTCCAC
TCACTGGGACTGTACCGCCAGGTGTCATGCGTCTCTGAGGTTTCTGATTAAAGGTTGTCTCGGTTT
CCCCA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_152718.2</u>
Summary:	May be a regulatory element in the beta-catenin signaling pathway and a target for chemoprevention of hepatocellular carcinoma.[UniProtKB/Swiss-Prot Function]
Locus ID:	220001
MW:	15.5