

Product datasheet for **SC205278**

SRPX (NM_006307) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SRPX (NM_006307) Human 3' UTR Clone
Symbol:	SRPX
Synonyms:	DRS; ETX1; HEL-S-83p; SRPX1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006307
Insert Size:	402 bp
Insert Sequence:	>SC205278 3'UTR clone of NM_006307 The sequence shown below is from the reference sequence of NM_006307. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GCCGAAATGAGCCAGACCTGTAACACCTGACATGATGGTTCCTCTCTTGGCAATTCCTCTTCATTGTCT
ACATAGTGACATGCACACGGGAAAGCCTTAAAAATATCCTTGATGTACAGATTTTATTTGTAATTTTAA
AAGTCTATTTTATTATGAGCTTTCTTTGCACTTAAAAATTAGCATGCTGCTTTTGTACTTGGAAGTGT
TTCAAAAAATTATATGACCATATTTACTCTTTCTAACTTTCTTTACTCCATCATGGCTGCTTGATTTTG
TAGAGAAATTAGAACCCATAACCATACACAGGCTATCAACATGTTATTCAATGTGACACCTAACTCTTT
TCTATTTTGTGTTTTTAAGTAAGACTTTTATTAATAAAACAAAATGTTTTGGAGCAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_006307.5

Summary: May be involved in phagocytosis during disk shedding, cell adhesion to cells other than the pigment epithelium or signal transduction.[UniProtKB/Swiss-Prot Function]

Locus ID: 8406

MW: 15.6